



TETRAGRAPH

PEDIATRIC **Sensor Placement**

Quantitative
Train-of-Four Monitor



Ulnar Nerve

Option 1

Grounding

Landmark

Hand


1 cm Below
the Wrist Crease

Foot

Medial Malleolus



Option 2



Posterior Tibial Nerve

Option 3



Stimulating Ulnar Nerve | Posterior Tibial Nerve



Visit [Senzime.com/support](https://www.senzime.com/support)

For more information, scan the QR code or
text our customer support line: (314) 788-7400.

*ADM Abductor Digiti Minimi muscle | AP Adductor Pollicis muscle | FHB Flexor Hallucis Brevis muscle

