

TetraGraph^{philips interface}

**Are you a
Philips user?**

Do you need the data
in your PDMS system?

**Get access to
EMG technology**

Data will be treated
as all Philips data



Data transfer
from TetraGraph
neuromuscular
monitor to
Philips IntelliVue
system

CONNECT.

VIEW.

TRANSFER.

CONNECT

Transfer NMT monitoring parameters; TOFR, TOFC, TOF bars, PTC, ST and Trend Graph

TetraGraph

TetraGraph Philips Interface

VIEW / TRANSFER

NMT data are displayed on the Philips monitor and will be integrated into existing infrastructure allowing upload of patient data to local PDMS (Patient Data Management System).

IntelliVue

EC10 Module



Learn more?
senzime.com

Order Information

SEN 2007 TetraGraph
Philips Interface

Supported Philips monitors

- Philips CMS monitor family. Software Release C or higher
- Philips IntelliVue MP40/50/60/70/80/90/MX400/450/500/550/600/700/750/800/850 Software Revision H15 or higher

Supported TetraGraph software versions

34b.22b.12a or higher

Compatible Philips IntelliBridge system

- PIIC iX IntelliVue Information Center iX (SW version B.0)
- PIC iX Patient Information Center iX (SW version C.0 or higher)

EC10 module

- OpenInterface driver version A.6 or higher.



SENZIME

senzime.com
info@senzime.com